

National Oceanography Centre, Southampton UNIVERSITY OF SOUTHAMPTON AND NATURAL ENVIRONMENT RESEARCH COUNCIL



# Carbonate Chemistry Facility for UK Ocean Acidification Research Programme

Eithne Tynan, Eric Achterberg, Cynthia Dumousseaud, Ziaul Chowdhury







- High quality measurements of dissolved inorganic carbon (DIC) and total alkalinity (TA).
  - Precision: ±0.1 % or better
  - Accuracy: CRMs (Prof. A. Dickson; Scripps)




#### Samples analysed up to March 2012:

Experimental samples	160
Oceanographic samples	700
Total	860



- High quality measurements of dissolved inorganic carbon (DIC) and total alkalinity (TA).
  - Precision: ±0.1 % or better
  - Accuracy: CRMs (Prof. A. Dickson; Scripps)







- High quality measurements of dissolved inorganic carbon (DIC) and total alkalinity (TA).
  - Precision: ±0.1 % or better
  - Accuracy: CRMs (Prof. A. Dickson; Scripps)





#### International inter-calibration exercise:

Involved in inter-laboratory comparison for sea water  $CO_2$  measurements (scheduled summer 2012)



- D366 UK cruise (summer 2011): Carbonate chemistry measurements for the bioassays





### Coming up:

JR271 Arctic cruise (summer 2012) JR274 SO cruise (winter 2012/13)

# Services of the facility

- Guide for best practices (including details on sampling, transport and use of the facility's online sample logging and receipt system)



Practical guide – UK Ocean Acidification Research Programme (UKOARP) Carbonate Chemistry Facility

Table of contents	
Introduction	4
1. Sampling and storage recommendations	
2. Total Dissolved Inorganic Carbon	
2.1. Coulometric titration	
2.2. Infrared detection	
3. Determination of Total Alkalinity	
3.1. VINDTA 3C	
3.2. Apollo SciTech AS-ALK2	
4. Precision and accuracy	
5. Calculation of the carbonate system parameters	
6. LIMS Software	
References	1

Practical guide available at http://www.noc.soton.ac.uk/ukoacf

For more information on recommended procedures:

<u>Dickson et al. 2007</u>: Guide to best practices for ocean  $CO_2$  measurements <u>EPOCA</u>: Guide to best practices for ocean acidification research and data reporting

## Services of the facility

- Training workshops at the National Oceanography Centre Southampton for sampling, storage, transport and data analysis



Workshop	Participants
8 <sup>th</sup> April 2011	20
12 <sup>th</sup> January 2012	7
2013	?


-Telephone and email technical support to the UK OA RP researchers

Prof. Eric Achterberg (eric@noc.soton.ac.uk; +44 23 8059 3199) Dr. Eithne Tynan (e.tynan@noc.soton.ac.uk; +44 23 8059 3641)

# Carbonate System Facility LIMS Software

#### Software available at: http://www.noc.soton.ac.uk/ukoacf

### Cost of analysis: £12 per sample (for both DIC and TA) for UK OA projects



Carbonate Facility LIMS website

# Carbonate System Facility LIMS Software

Username and password protected

#### Input by user:

- sample information (ID, number...)
- ancillary data: salinity, temperature, pressure, phosphate and silicate concentration
- optional: pH (including pH scale), pCO<sub>2</sub>

### **Output by Carbonate System Facility (once samples analysed):**

- analysed TA and DIC
- calculated pH (total scale)
- calculated pCO<sub>2</sub>
- carbonate and bicarbonate ion concentration
- calcite and aragonite saturation state

# Carbonate System Facility LIMS Software

Username and password protected

#### Input by user:

- sample information (ID, number...)
- ancillary data: salinity, temperature, pressure concentration
- optional: pH (including pH scale)

### Output by Carbonate S

- analysed TA and
- calculated
- calculated
- carbonate ar carbonate ion concentration
- calcite and aragonite saturation state

d silicate

once samples analysed):